

Spa Range

UBS1000 Urinal Bottle - Male urinal with handle, hospital style

Long Life Re-usable Plastic Male Urinal Bottle for use in Hospitals, Nursing, Rehabilitation, Elderly Care, and Community Care



The Spa Range is a range of re-usable living aids designed for use in Hospital Wards, Nursing and Care Homes, Rehabilitation and Personal Use at home.

Our Spa Range is an attractive, stylish and functional group of products which includes drinking cups, beakers, jugs, vomit bowls, medicine cups, urinals, bedpans and commodes.



Product Benefits

- Lightweight, strong and durable.
- Secure fitting cap to prevent leaks and spills and reduces unpleasant odours and risk of cross infection.
- Graduated with clear, easy to read measurements.
- Easy grip handle - particularly suitable for those with limited hand movement.
- Hollow handle for improved cleaning results.
- Can be cleaned in automatic washer disinfectors or hand washed.
- Smooth material surfaces ensure ease of use and ease of cleaning.
- Guaranteed for a minimum of 1000 wash cycles.
- Dual graduated scale, allowing contents to be measured when bottle is placed vertically or horizontally.

Product Specification

- Male Urinal Bottle - Graduated, 1200ml
- Pack Size 45
- Pack Weight 7.5kg
- The following standards are applicable:
 - **ISO 9001:2008**
Quality Management System
 - **ISO 13485:2003**
Medical Devices - Quality Management System
 - **BS3215:1995**
Re-usable Portable Urinals
- Manufactured from Polyethylene (HDPE)

Guide to Using Spa Products from Warwick SASCo

Usage

- For male use only.
- Can be used in seated, standing or supine positions.
- Close cap immediately after use to prevent spillage and odours.
- Empty after each use.

Cleaning

- **For Hospital Use** - Wash before first use. Clean in automatic washer disinfectors compliant to ISO 15883. Maximum water temperature 90°C.
- For further processing information, refer to machine manufacturers instructions and local hospital policy.
- **For Personal Use** - Wash before first use in hot water with detergent.
- Empty contents immediately after use. Wash in hot water with detergent after each use and allow to dry. .

Storage

- This product has a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.