


Comfy Gel



Comfy Gel

- EZ Feel Gel Covered by a Sheepskin Cover
- Cover is made from Medical Sheepskin
- Has the Comfort Advantage of Gel Cushion & Sheepskin
- Available in 1 inch (25 mm) or 3 Inch (75 mm) thick
- 46 cm x 40 cm
- Code: (25 mm Thick) C8001
(75 mm Thick) C8002


EZ Feel Gel



EZ Feel Gel

- High Quality Polyurethane Gel
- Excellent Pressure Distribution & Comfort
- Gel will not Migrate, Leak, Evaporate, Harden or Dry out if Punctured
- Resistant to Bacteria and Bodily Fluids
- Easily Cleaned with Soapy Water or Disinfectant
- 46 cm x 40 cm
- Water Proof
- Code: NC1816


EZ Feel Gel - Lite 3



EZ Feel Gel - Lite 3

- EZ Feel Cushion on a High Density Polyurethane Foam
- 75mm Foam Under Support
- Assists in Preventing Pressure Sores
- Provides Excellent Immersion Stability and Pressure Distribution
- Suitable for Clients with Moderate to High Risk of Skin Breakdown
- 46 cm x 40 cm
- Can be combined with Sheepskin Cover
- Code: EZ1816-3


Ultra Gel / Air



Ultra Gel / Air

- Gel/Air Design reduces Overall Weight of Cushion
- Gel/Air provides a Gel Characteristic Combined with Air
- Quilted Gel Provides Outstanding Pressure Relief
- 46 cm x 40 cm
- Can be combined with Sheepskin Cover
- Code: NAL1816

Ultra Gel / Air - Lite 3



Ultra Gel / Air - Lite 3

- Gel/Air Design Reduces Overall Weight of Cushion
- Gel/Air provides a Gel Characteristic Combined with Air
- Quilted Gel Provides Outstanding Pressure Relief
- 75mm Foam Moulded Under Support
- 46 cm x 40 cm
- Can be combined with Sheepskin Cover
- Code: UL1816-3 (46cm x 40cm) or UL2018-3 (50cm x 46 cm)



Stream Lite - PU



Stream Lite - PU

The Stream Line - PU cushion produces skin protection with the positioning and stability for more developed function. This contoured shape positions the buttocks in the wheelchair, and higher cells in the middle of the cushion supplying better stability.

- Code: 1R1816-SL

Cushion Air - PU



Cushion Air - PU

The Cushion Air - PU delivers a high level of skin protection and stability while supplying the proper amount of postural support.

2 Way stretchable cover enables the cushion cells to conform to the individual's body shape without diminishing product effectiveness or therapeutic value.

CushionAir - PU consists of dual compartments and can be extended to 4 compartments if required, which react independent to inflation and pressure.

- Code: 1 Chamber: 1R1816-4
2 Chamber: 2R1816-2

Standard Cover



Standard Cover (supplied with all cushions)

- Stretchable, low-shear top enables the individual cells to conform to different body shapes without diminishing the cushion's therapeutic value.
- Non-slip bottom helps secure the cushion to the chair.
- Four-sided fishnet design increases ventilation. Side eyelets provide easy access to valves for inflation.

Sheepskin Cover - Optional



Sheepskin Cover - optional

- Washable Mediwool
- Available in 25mm thick and 75mm thick

Code: 25mm: CUSH1
75mm: CUSH2

